



Art.-No. 62390

Mouse Trap with bait, Set of 2

Rodent monitoring



Pull off the protective cap, set the trap by a wall, tension it and it is ready. As this plastic trap is pre-baited, it is ready for immediate use. To remove a caught rodent, simply press the tensioning bar again. Touching the mouse is not necessary. As the trap is made of plastic, it can be rinsed for hygienic re-use.

We recommend placing mouse traps by the wall. Mice prefer to run there, as they feel safer than in the open.

N-98146

Use biocides safely. Always read the label and product information before use.



Additional information

- **Dimensions:** 4.8 x 5.5 x 10.5 cm (W x H x D)
- **Packaging dimensions:** 12.9 x 17 cm (W x H)
- **Manuals included:** DE, GB, FR, ES, IT, NL
- **Product images, scenery images**
- **Packaging unit:** 100 units per package

Easy and quick to set up!

